



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Double Hook S.S Saline Stand

**Product Code :**  
SAMCOMEDICAL120001A013



**Description :**

Double Hook S.S Saline Stand with IV Stand Economy

**Technical Specification :**

Saline Stand/IV Stand - Economy M.S. tubular pipe. Five legged heavy polypropylene base mounted on 5 cms castors. Double hook S.S. adjustable IV rod with a bakelite knob. Pre-treated and epoxy powder coated. Knockdown construction



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),