



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Stretcher On ERW Tubular Trolley

Product Code :
SAMCOMEDICAL120001A017



Description :

Stretcher On ERW Tubular Trolley

Technical Specification :

Stretcher On Trolley Size : 210L x 56W x 82H cms. Frame work made of ERW tubular pipe of 31 mm. dia mounted on four 15 cms. dia. castors (2 with brakes). Removable stretcher top made of S.S. Sheet pressed on tubular frame. Brackets on alternate ends for IV drip rod. Knockdown construction



Samco Medical

Website: www.samcomedical.com,