



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Portable Battery Operated Suction Unit (Eurovac B)

**Product Code :**  
SAMCOMEDICAL120001A023



**Description :**

Portable Battery Operated Suction Unit (Eurovac B) Capacity: -600 mm Hg to 10 Vacuum, 30 ltrs/ min. Power: Operates on 220 V AC and in-built maintenance free 12 V, 7.5 Amp rechargeable battery(110 V on request)  
Vacuum Gauge: 6.25 cms. 0-760mm Hg calibration. Jars: Wide mouthed autoclavable 2 ltr. collection jar (Polycarbonate), safety jar with mechanical/electronic overflow motor shut-off system and autoclavable bacterial filter. Pump Type: Oil free diaphragm pump Housing: Engineering moulded ABS Plastic Cabinet with electronic overflow protection safety jar Dimensions & Weight: 23 X 42 X 36cms & 11Kg. Battery Backup: 120 minutes on full charge. In built maintenance free battery. Tubing: 8mm ID x 2mtr. (Silicone)



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),