



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Samcomedical Suction Unit In Stand AC with DC

Product Code :
SAMCOMEDICAL120001A024



Description :

Samcomedical Suction Unit In Stand AC with DC

Technical Specification :

Suction Unit - INSTAD AC/DC Housing: ABS Plastic Moulded & 3 Stage LED Indicator for Battery Status
Capacity: -600 mm Hg 10 , 21 Ltrs/ min Pump Type: Piston Pump Jars: 1 x 1 Ltr Polycarbonate Jar Filter:
Bacterial Filter Autoclavable/Reusable Tubing: 8 mm ID x 2 mtr. (Silicon) Vacuum Gauge: 5 cm dia 0-760 mm
Hg calibration. Power: 220 V AC, 50 / 60 Hz, 70Watt , 12 V DC (Internal Rechargeable Battery 12V 4.5Amp)
Battery: Lead Acid (Ni-cd Optional) Noise Level: 55 dB A 3 Dimension & Weight: 39 x 22 x 22 cms, 6 Kg.
Overflow Safety: Mechanical Accessories: Car Charger & Power Adaptor



Samco Medical

Website: www.samcomedical.com,