



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Oxygen Concentrators and Cylinder

**Product Code :**  
SAMCOMEDICAL8001A047

**Description :**

Oxygen Concentrators and Cylinder

**Technical Specification :**

Oxygen Concentrators and Cylinder Oxygen concentration: ?90% Oxygen flow: 1~5L/min Output pressure: 30-70kPa Overrunning pressure alarm Noise: ?52dB(A) Power supply: AC220ï¿½22V, 50ï¿½1Hz Input power: ?500VA



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),