



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Oxygen Concentrators Nebulizing Installation

Product Code :
SAMCOMEDICAL8001A050

Description :

Oxygen Concentrators Nebulizing Installation

Technical Specification :

Oxygen Concentrators Nebulizing Installation It apply zeolite molecular sieve as shift, adopt advanced physical theory(PSA) to separate oxygen and nitrogen in air. At the same time filter the bad substance in air, so as to get pure oxygen for medical use Oxygen Concentration : 93%±3%. Oxygen flux: 1~3L/min. Output pressure: 0.04~0.05Mpa. Overrunning pressure alarm. Operation noise: 49dbA. Power supply: AC220V±20V, 50Hz±1 Hz. Input power: 350VA



Samco Medical

Website: www.samcomedical.com,