



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
C-Arm Compatible, Fully Electromatic O.T. Table

Product Code :
SAMCOMEDICAL9007A004



Description :

C-Arm Compatible, Fully Electromatic O.T. Table Technical Constraints Top Dimension : L 2032 xW 550mm
Height Adjustment : 787mm- 1037mm Table Top Sliding : 250mm(optional) Trendelenburg / Reverse : 30° / 25° Lateral Tilt : 20° / 20° Kidney Elevator : 150mm Back Rest (up / down) : 80° / 25° Leg Rest (up / down) : 15° / 90° Head Rest (up / down) : 20° / 60° Power Supply : 24V DC Battery Backup (Optional at extra price) : 2 - 3 hours Standard Accessories : Anaesthetic screen : 1 pc. Shoulder supports : 2 pcs. Lateral supports : 2 pcs. Arm rest (S.S.top) : 1 pc. Knee crutches : 2 pcs Overview : This Electromatic O.T Table is suitable for all kind of surgeries, i.e. general surgery, vascular surgery, cardio, neurology, urology, gynecology, proctology, laparoscopy, traumatic surgery, dental surgery, plastic surgery and others Key Structure: The five section table top is X-ray translucent and is equipped with removable, antistatic 50mmthick mattresses. The base & column of the table is made of high quality medical stainless steel (Grade-304). Non-reflecting surface, antibacterial & is easy to clean. The table can provide with a top sliding of 250 mm, Dual Leg Section & Battery backup as per consumer requirement. (Optional) Height Adjustment control by foot controller. Head and Leg sections are detachable & interchangeable. Longitudinal Sliding Top, Lateral Tilts, Trendelenburg, Flex-Reflex, Chair Position and are realized by Remote control



Samco Medical

Website: www.samcomedical.com,