



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Electro Mechanical Operating Table Remote Controlled Model

Product Code :
SAMCOMEDICAL9007A006



Description :

Electro Mechanical Operating Table Remote Controlled Model Specifications/Details Four section eccentrically positioned radio-translucent table top suitable for C-arm image intensifier. Operating positions : Height adjustment, lateral tilt, trendelenburg, reverse trendelenburg and back section are precisely and smoothly controlled by hand switch through imported electro mechanical actuators. Head & foot section manually operated by the means of ratchet system. Complete with standard accessories and mattress. The electro-mechanical driving mechanism assures smooth and quiet movements of the table-top. Operating positions including height adjustment, trendelenburg, reverse trendelenburg, lateral tilt and back section can be controlled quickly and accurately with the help of lightweight hand control unit. The excentrically positioned table top, without cross members, to optimize image intensifier accessibility for both A/P & lateral exposures. The table top is made of radio-translucent material for fluoroscopic capability. Interchangeable head & leg sections and manually operated by the means of precise ratchet system. Base and column covers are made of stainless steel for easy cleaning and hygiene Standard Accessories L-shaped anaesthetic frame : 1 pc. Shoulder support with pad : 1 pair Arm board with S.S. top : 2 pcs. Lateral support with pad : 1 pair. Wrist strap : 1 pair. Knee crutches goepel type : 1 pair. Water proof rubber mattress : 1 Set



Samco Medical

Website: www.samcomedical.com,

