

Product Name :
Hydraulic Surgical Operating Table With Stainless Steel
Embedded Base And All Accessories

Product Code :
SAMCOMEDICAL9007A007



Description :

Hydraulic Surgical Operating Table With Stainless Steel Embedded Base And All Accessories
Specifications/Details Five section stainless steel top. Operating positions : Back section, trendelenburg, reverse trendelenburg, and lateral tilt etc., are precisely and smoothly controlled from either side of the table by working of different handles. With stainless steel embedded base and all accessories of stainless steel. Complete with standard accessories and mattress. Precision made hydraulic system with a single lever action hydraulic lift controlled by a foot pedal. Five section stainless steel top gives anatomic adaptation to the patient body. Built-in Kidney bridge and large perineal cut-out for drainage. All the operating positions including proctology, lithotomy, chair, flex, reflex, trendelenburg and anti trendelenburg, lateral tilt etc. Can be obtained from either side of the table by working of different handles. Precise and easy to operate gear systems for smooth and accurate table movements. Lithotomy positioning : By simply removing the leg-end section and adding stirrups (optional) and drainage tray (optional), the table easily converts from supine or prone position to lithotomy position. Heavy cast iron steel base provides stability to the entire table and has toe space and foot rails for comfortable footrest for the surgeons. Mounted on four heavy duty swivel castors with floor locking device. Fitted with side rails and adjustable clamps to fix the accessories at any desired position
Standard Accessories L-shaped anaesthetic frame : 1 pc. Shoulder support with pad : 1 pair. Intravenous arm board with s.s. Top : 1 pc. Padded side supports : 1 pair. Wristlets : 1 pair. Lithotomy leg holder with straps : 1 pair. Foot rest with ss top : 1 pc

Website: www.samcomedical.com,

