



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Gearmatic Hydraulic OT Table With Five Section
Translucent Top.

Product Code :
SAMCOMEDICAL9007A008



Description :

Gearmatic Hydraulic OT Table With Five Section Translucent Top. The top segments are translucent for X-rays and enable taking pictures and monitoring a patient with a C-arm. The positions in five section table top are adjusted manually by mechanical drive. It's safe and robust mechanism prevents breakdown during an operation# Key Structure: With regular human strength, up-down table position can be fixed by stepping foot pedal downward. The base is eccentric for trouble-free C-ARM movements. Leg & Head section are interchangeable. Base & Column covered with stainless steel Technical Constraints Top Dimension : L 1905 xW 553mm Height Adjustment : 787mm- 1042mm Trendelenburg / Reverse : 30° / 25° Lateral Tilt : 20° / 20° Kidney Elevator : 150mm Back Rest (up / down) : 80° / 25° Leg Rest (up / down) : 15° / 90° Head Rest (up / down) : 20° / 60° Standard Accessories : Anaesthetic screen : 1 pc. Shoulder supports : 2 pcs. Lateral supports : 2 pcs. Arm rest (S.S.top) : 1 pc. Knee crutches : 2 pcs



Samco Medical

Website: www.samcomedical.com,