

Product Name :
Side End Control Super Deluxe Hydraulic Operating Table

Product Code :
SAMCOMEDICAL9007A009



Description :

Side End Control Super Deluxe Hydraulic Operating Table Versatile Super Deluxe O.T. Table is superior for general surgical procedures. Accurate positions can be achieved with easy mechanical operation. Five section stainless steel top & heavy base and column enhances the durability and hygiene Key Structure: The five section Stainless Steel table top is equipped with removable, antistatic mattresses. All the positions are realized by a mechanical driver. The table is furnished with hard, robust stainless steel shell which provides the sense of dirt free & sanitized operations Technical Constraints Top Dimension : L 1905 xW 553mm Height Adjustment : 787mm- 1042mm Trendelenburg / Reverse : 30° / 25° Lateral Tilt : 20° / 20° Kidney Elevator : 150mm Back Rest (up / down) : 80° / 25° Leg Rest (up / down) : 15° / 90° Head Rest (up / down) : 20° / 60° Standard Accessories : Anaesthetic screen : 1 pc. Shoulder supports : 2 pcs. Lateral supports : 2 pcs. Arm rest (S.S.top) : 1 pc. Knee crutches : 2 pcs



Samco Medical

Website: www.samcomedical.com,