

Product Name :
C-Arm Image Intensifier Hydraulic OT Table

Product Code :
SAMCOMEDICAL9007A010



Description :

C-Arm Image Intensifier Hydraulic OT Table Key Structure The 5 section table top contains radiographic sheet of 12 mm and is equipped with assembled polyurethane 50mmthick mattresses. Easy to disassemble & cleanse. Strong & rust free base of stainless steel is provided to make the working conditions pleasant for the doctors & nurses. Head / Leg section are interchangeable. Auto floor locking device Technical Constraints Top Dimension : L 1890 xW 510mm Height Adjustment : 760mm- 990mm Trendelenburg / Reverse : 30° / 30° Lateral Tilt : 20° / 20° Kidney Elevator : 150mm Back Rest (up / down) : 80° / 25° Leg Rest (up / down) : 15° / 90° Head Rest (up / down) : 20° / 90° Standard Accessories Anaesthetic screen : 1 pc. Shoulder supports : 2 pcs. Lateral supports : 2 pcs. Arm rest (S.S.top) : 1 pc. Knee crutches : 2 pcs Key Structure : Provide 360-degree rotation positioning of the head system. Adopts all major headrest and skull clamps to the Mayfield three head pin base unit. Exclusive spring loaded mechanism provides assured stability and a secured lock. The skull clamp has the additional feature of a swiveling c-arm. This swiveling capacity greatly facilitates positioning the pins on the head around nerves, vessels, and thin bone in the temporal region. For better result, the arm is locked by control knob and the width is adjusted by ratchet mechanism. The three cone-shaped pins available in adult and child sizes are interchangeable and may be removed for sterilization Key Structure : This attachment, originally designed for transsphenoidal hypophysectomy is used for all prone position procedures. It is attached to the adjustable basic unit by means of transitional member and the gardner clamp adaptor. It is padded with soft, pliable sponge rubber for equalized weight distribution and incorporates the new quick -change pads. The horseshoe width can be quickly adjusted to accommodate various head sizes. A removable, adjustable pulley bar is included if traction is desired

Samco Medical

Website: www.samcomedical.com,

