



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Electric Operating Table - Orange Series

Product Code :
SAMCOMEDICAL9007A013



Description :

Electric Operating Table - Orange Series Technical specifications Panel material : Carbon fiber Table Top (Length X Width) : 2010 $\frac{1}{2}$ 10 mmX535 $\frac{1}{2}$ 10mm Height Adjustment : 995 $\frac{1}{2}$ 20 mm-695 $\frac{1}{2}$ 20mm Height range : 300 $\frac{1}{2}$ 20 mm(electric) Reversed Trendelenburg : 28 $\frac{1}{2}$ 2 $\frac{1}{2}$ (electric) Trendelenburg : 28 $\frac{1}{2}$ 2 $\frac{1}{2}$ (electric) Lateral Tilt : 23 $\frac{1}{2}$ 2 $\frac{1}{2}$ (electric) Head Section up/down : 45 $\frac{1}{2}$ (manual)/ 90 $\frac{1}{2}$ (manual) Back Section up/down : 60 $\frac{1}{2}$ (manual)/ 35 $\frac{1}{2}$ (manual) Leg Section up/down : 35 $\frac{1}{2}$ (manual)/ 90 $\frac{1}{2}$ (manual) Translation distance : 300mm (manual) Bulid-in battery : 24V 15AH Kidney Section Elevation : 120mm Supply : AC 220 $\frac{1}{2}$ 10%,50HZ



Samco Medical

Website: www.samcomedical.com,