



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Electrically Operated OT Table with Double Layer Table
Top Suitable for X-Ray Examintaion

Product Code :
SAMCOMEDICAL9007A016



Description :

Electrically Operated OT Table with Double Layer Table Top Suitable for X-Ray Examintaion Standard Features: Suitable for surgeries of head, neck, pleuroperitoneal cavity, perineum, limbs, gynecology and obstetrics department, urology thoracic surgery etc. Double layer table top suitable for x-ray examination. Framework & accessories made of stainless steel. Mattress is resistant against acids, alkaline solutions and is anti-bacterial & can be cleaned easily. Various positions like height adjustment, trendelenburg, reverse trendelenburg, lateral tilt, back rest all controlled by linear power device. Equipped with auto chargeable device & battery back up. Leg sections are detachable. In leg section, leg plates are dismountable and separable Technical specifications Table Top (Length x Width) : 2050mm x 480mm Height Adjustable : 700mm-950mm Trendelenburg : $\pm 15^{\circ}$ / $\pm 25^{\circ}$ Reverse Trendelenburg : $\pm 15^{\circ}$ / $\pm 25^{\circ}$ Lateral Tilt : $\pm 15^{\circ}$ / $\pm 20^{\circ}$ / $\pm 30^{\circ}$ Back Rest (up/down) : $\pm 15^{\circ}$ / $\pm 80^{\circ}$ / $\pm 15^{\circ}$ Head Rest (up/down) : $\pm 15^{\circ}$ / $\pm 80^{\circ}$ / $\pm 80^{\circ}$ Kidney Bridge Elevation : ± 120 mm Leg Rest adjustment : $\pm 90^{\circ}$ / $\pm 15^{\circ}$ Power : 220V / 50HZ Foot rest can be kept horizontally or perpendicularly to the Leg Section. Standard Accessories : Knee Crutch : 2 pcs Body Strap : 2 pcs Arm Rest : 2 pcs Anesthesia Rack : 1 pc. Controller : 1 pc



Samco Medical

Website: www.samcomedical.com,