

Product Name :
Electrically Operated OT Table

Product Code :
SAMCOMEDICAL9007A017



Description :

Electrically Operated OT Table Standard Features: Suitable for surgeries of urinary, obstetrical, gynaecological and so on. It is an ideal universal operating table. Double layer table top (made of paper plywood) suitable for x-ray examination. It is quite convenient for various x-ray positions. Mattress base frame is made of aluminium. Accessories, base shell and column shell all made of stainless steel. Mattress is resistant against acids, alkaline solutions and is anti-bacterial & can be cleaned easily. Hydraulic pump and electro-magnetic valve are imported and of very good quality. Leg section is detachable. In leg section, leg plates are dismountable and separable and can be folded down 90°. It is suitable for urinary, obstetrics and Gynaecology surgery

Technical specifications

Table Top (Length x Width) : 2050mmx480mm Height Adjustable : 700mm-950mm Trendelenburg : 15°/20° Reverse Trendelenburg : 15°/20° Lateral Tilt : 15°/18° / 18° Back Rest (up/down) : 15°/75° / 20° Head Rest (up/down) : 15°/75° / 90° Kidney Bridge Elevation : 120mm Leg Rest adjustment : 90° Power : 220V / 50HZ Foot rest can be kept horizontally or perpendicularly to the Leg Section. Standard Accessories : Knee Crutch : 2 pcs Body Strap : 2 pcs Arm Rest : 2 pcs Anesthesia Rack : 1 pc. Controller : 1 pcs



Samco Medical

Website: www.samcomedical.com,