



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Head End Controlled Universal OT Table

**Product Code :**  
SAMCOMEDICAL9007A018



**Description :**

Head End Controlled Universal OT Table Standard Features: Suitable for thoracic surgery, abdominal surgery, ophthalmology, otorhinolaryngology. Accessories and framework of the table made of stainless steel. Height is adjusted with foot operated oil pump. Other various movements were handled by mechanical transmission head. Mattress is resistant against acids, alkaline solutions and is anti-bacterial & can be cleaned easily

Technical specifications Table Top (Length x Width) : 2100mm x 480mm Height Adjustable : 800mm-1045mm  
Trendelenburg :  $\pm 15^\circ$  Reverse Trendelenburg :  $\pm 15^\circ$  Lateral Tilt :  $\pm 25^\circ$  /  $25^\circ$  Back Rest (up/down) :  $\pm 85^\circ$  /  $10^\circ$  Head Rest (up/down) :  $\pm 45^\circ$  /  $50^\circ$  Kidney Bridge Elevation : It can be folded upto  $90^\circ$  Leg Rest adjustment (down) :  $90^\circ$  Standard Accessories : 2 knee crutches, 2 body strap, 2 armrests & 1 anaesthesia rack



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),