



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Luxury Multifunctional Operating Table

**Product Code :**  
SAMCOMEDICAL9007A020



**Description :**

Luxury Multifunctional Operating Table Standard Features: Suitable for surgeries of chest, abdomen, gynaecology and obstetrics, orthopaedics. Accessories and framework of the table made of stainless steel. Height is adjusted with foot operated oil pump. Back section and leg section are operated by imported gas spring technology. Mattress is resistant against acids, alkaline solutions and is anti-bacterial & can be cleaned easily

Technical specifications Table Top (Length x Width) : 2000mm x 480mm Height Adjustable : 650mm-890mm  
Trendelenburg :  $\pm 28^{\circ}$  Reverse Trendelenburg :  $\pm 28^{\circ}$  Lateral Tilt :  $\pm 20^{\circ}$  /  $20^{\circ}$  Back Rest Adjustable :  $\pm 75^{\circ}$  /  $10^{\circ}$  Head Rest (up/down) :  $\pm 90^{\circ}$  /  $90^{\circ}$  Kidney Bridge Elevation : 100mm Leg Rest adjustment (down) :  $90^{\circ}$  (Detachable) Standard Accessories: Knee Crutch : 2 pcs Body Strap : 2 pcs Arm Rest : 2 pcs Anesthesia Rack : 1 pc



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),