



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Hydraulically Operated Multi-Purpose OT Table

Product Code :
SAMCOMEDICAL9007A021



Description :

Hydraulically Operated Multi-Purpose OT Table Standard Features: Suitable for thoracic, abdominal, ophthalmology, otorhinolaryngology, obstetrics and gynecology, urinary surgery and orthopedic surgery. Lifting of this operating table is controlled by foot oil pump and rest of the movements are controlled by mechanical transmission side manipulation Technical specifications Table Top : 2100mm x 480mm Height Adjustment : 1045mm-800mm Reversed Trendelenburg : $\pm 25^{\circ} \frac{1}{2}$ Trendelenburg : $\pm 25^{\circ} \frac{1}{2}$ Lateral Tilt : $\pm 20^{\circ} \frac{1}{2}$ Head Rest down : $\pm 20^{\circ} \frac{1}{2}$ Head Rest size : 275 x 310mm Vertical of oil pump : 240mm Back Section size : 500 x 480mm Leg Section size : 500 x 480mm Seat section size : 450 x 480mm



Samco Medical

Website: www.samcomedical.com,