

Product Name :
Electro Hydraulic Operated Multi Purpose OT Table

Product Code :
SAMCOMEDICAL9007A022



Description :

Electro Hydraulic Operated Multi Purpose OT Table Standard Features: Suitable for thoracic surgery, abdominal surgery, orthopaedics, cerebral surgery, ophthalmology, otorhinolaryngology, obstetrics and gynecology and urinary surgery. X-ray permeable table top cassette guide rail available under the table top. Framework & accessories made of stainless steel. Mattress is resistant against acids, alkaline solutions and is anti-bacterial & can be cleaned easily. Various positions like height adjustment, trendelenburg, reverse trendelenburg, lateral tilt, back rest all controlled by electric hydropneumatic device. Head & Leg sections are detachable. In leg section, separable & dismountable leg plates are optional (at extra cost) Technical specifications Table Top (Length x Width) : 1960mm x 480mm Height Adjustable : 690mm-930mm Trendelenburg : $\pm 25^{\circ}$ Reverse Trendelenburg : $\pm 20^{\circ}$ Lateral Tilt : $\pm 18^{\circ}$ / 18° Back Rest (up/down) : 75° / 10° Head Rest (up/down) : 45° / 50° Kidney Bridge Elevation : 780mm Leg Rest adjustment : 90° Detachable Standard Accessories : Knee Crutch : 2 pcs Body Strap : 2 pcs Arm Rest : 2 pcs Anesthesia Rack : 1 pc. Controller : 1 pcs



Samco Medical

Website: www.samcomedical.com,