



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Electric Hydraulic Operation Table

**Product Code :**  
SAMCOMEDICAL9007A023



**Description :**

Electric Hydraulic Operation Table Standard Features: Suitable for thoracic surgery, abdominal surgery, orthopaedics, cerebral surgery, ophthalmology, otorhinolaryngology, obstetrics and gynecology and urinary surgery. X-ray permeable table top cassette guide rail available under the table top. Framework & accessories made of stainless steel. Integrated with rechargeable battery. "One Key reset" makes it convenient and easy to position the table after operation

Technical specifications

Table Top (Length x Width)	: 2020mm x 520mm
Height Adjustable	: 1100mm-700mm
Trendelenburg	: ?30°½
Reverse Trendelenburg	: ?30°½
Lateral Tilt	: ?20°½
Head Section up/ down	: ?40°½ / 90°½
Back Section up/ down	: ?80°½ / 45°½
Leg Section down	: ?90°½
Leg Section Opened	: ?90°½
Kidney Section Elevation	: 128mm
Build-in battery	: AC 220V½10%,50HZ



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),