



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Microprocessor Controlled OT Table With C Arm
Compatible Top

Product Code :
SAMCOMEDICAL9007A024



Description :

Microprocessor Controlled OT Table With C Arm Compatible Top Standard Features: The table is constructed with 5 movable parts, and suitable for all kinds of surgeries as well as special operative positions. C-arm compatible table top. Its unique design and large space under the table gives adequate space to install C arm X-Ray apparatus and other device. Equipped with X-ray rack and cassette specifically for C arm users. Equipped with microprocessor based controller and safety lock key protecting from mishandling. Set one-key restoration device. Structure made of stainless steel. Mattress is resistant against acids, alkaline solutions and is anti-bacterial & can be cleaned easily. Head & Leg sections are detachable. In leg section, leg plates are dismountable and separable. Generous, nice looking design is such that it give surgeons stable gravity centre and decreases fatigue & stress during operation Technical specifications: Table Top (Length x Width) : 2000mm x 520mm Height Adjustable : 650mm-950mm Trendelenburg : $\pm 15^{\circ}$ / $\pm 28^{\circ}$ Reverse Trendelenburg : $\pm 15^{\circ}$ / $\pm 28^{\circ}$ Lateral Tilt : $\pm 15^{\circ}$ / $\pm 18^{\circ}$ Back Rest (up/down) : $\pm 15^{\circ}$ / $\pm 75^{\circ}$ / $\pm 40^{\circ}$ Head Rest (up/down) : $\pm 15^{\circ}$ / $\pm 80^{\circ}$ / $\pm 80^{\circ}$ Leg Rest : downward 90° / $\pm 15^{\circ}$, Expanding 180° / $\pm 15^{\circ}$ Power : 220V / 50HZ Standard Accessories: Knee Crutch : 2 pcs Body Strap : 2 pcs Arm Rest : 2 pcs Anesthesia Rack : 1 pc. Controller : 2 pcs



Samco Medical

Website: www.samcomedical.com,