



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Taylors Brace Dorso Lumbar Brace

**Product Code :**  
NMBSB10003

**Description :**

Taylors Brace Dorso Lumbar Brace

**Technical Specification :**

Taylors Brace/Dorso Lumbar Brace ize S, M, L, XL Feature: Elastic abdominal panels provide rigid compression Thoracic straps provide comfortable fixation of brace to the body. Adjustable para-spinal and lateral bars are precontoured to keep spine in it's right position Clavicle straps with cotton stockinet for right shoulder posture Foam pads for good cushioning Application: Used for post-operative cases Mild thoraces and Lumber injuries Herniated Disc problems



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),