



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Cervical Collar Hard Adjustable

Product Code :
SAMCOMEDICAL100001A008

Description :

Cervical Collar Hard Adjustable

Technical Specification :

Cervical Collar Hard Adjustable Size : S, M, L Feature : High-density polythene sheet Proper foam padding for comfort Double layer facilitates height adjustment Holes for ventilation Velcro closure for easy application
Application : For cervical spine immobilization Used during Pre and post surgical cases Provides rigid support and strong immobilization



Samco Medical

Website: www.samcomedical.com,