



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Philadelphia Collar

**Product Code :**  
SAMCOMEDICAL100001A010

**Description :**

Philadelphia Collar

**Technical Specification :**

Philadelphia Collar Size: S, M, L, XL Feature: Two Piece design with proper velcro provision for comfortable fitting. Holes for ventilation Provides Preformed Chin support Front opening for tracheotomy Washable foam to maintain hygiene conditions Application: Pre and Post operative cases Severe Spondylosis Rigid support to the mild and lower cervical spine



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),