



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Knee Support with Hinges

**Product Code :**  
SAMCOMEDICAL100005A003

**Description :**

Knee Support with Hinges

**Technical Specification :**

Knee Support with Hinges Size : S, M, L, XL Feature : Wide elastic base with velcro / buckle combination for better compression Patellar opening provided to avoid pressure at patella Two lateral hinges provided for rigid compression Application : Rigid compression support for strains and sprains Post operative support Can be used for arthritis patients. Sports injuries



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),