



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Tubular Knee Support with Center Hole

**Product Code :**  
SAMCOMEDICAL100005A004

**Description :**

Tubular Knee Support with Center Hole

**Technical Specification :**

Tubular Knee Support with Center Hole Size : S, M, L, XL Feature : Soft yarn with elastic combination for mild compression Heat resistant rubber used Patellar opening provided to avoid pressure at patella Application : Mild support to injured knee in cases of strains and sprains Soft tissue injuries For stiff swollen and painful knee conditions



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),