



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Range of Motion Knee Brace

**Product Code :**  
SAMCOMEDICAL100005A012

**Description :**

Range of Motion Knee Brace

**Technical Specification :**

Range of Motion Knee Brace Size : Universal Feature : Dial screw to give desired angle for range of motion  
Foam padded splints to give comfortable fitting Application : Pre and post operative and rehabilitation care  
Anterior cruciate ligament (ACL) and Posterior cruciate ligament tear (PCL) Used for Total Knee Replacement (TKR) rehabilitation care



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),