



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :

Sphygmomanometer Aneroid Round Mobile

Product Code :

SAMCOMEDICAL6001A015

Description :

Sphygmomanometer Aneroid Round Mobile

Technical Specification :

Sphygmomanometer Aneroid Round Shaped Mobile, ABS Round dial face adjustable for easy reading, stand on aluminum tube with 5 plastic castors, metal basket, standard latex bulb, inflation valve with spring, standard end valve, PVC tubing. Sphygmomanometer Aneroid Round Mobile, Diagnostic Equipments Exporters, Diagnostic Equipments Manufacturers



Samco Medical

Website: www.samcomedical.com,