



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Cervical Traction Kit (Sitting)

**Product Code :**  
SAMCOMEDICAL100006A002

**Description :**

Cervical Traction Kit (Sitting)

**Technical Specification :**

Cervical Traction Kit (Sitting) Size : Universal Feature : Padding on head halter to give comfortable wear  
Spreader bar and pulley system provision Heavy-duty weight bag with separate graduations for water and sand  
Heavy-duty strings provided in the kit. Application : Cervical spine traction in sitting position For short time  
cervical traction in muscle spasm Cervical spine dislocation for continuous



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),