



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Skin Traction Kit

Product Code :
SAMCOMEDICAL100006A003

Description :

Skin Traction Kit

Technical Specification :

Skin Traction Kit Size : Child, Adult Feature : Foam padded heavy-duty fabric and crepe bandage provided to wrap around for below knee skin traction Heavy-duty strings provided Availability of Pulley system Heavy-duty weight bag with separate graduations for water and sand. Application : Below knee traction to relieve pain in the back and knee region



Samco Medical

Website: www.samcomedical.com,