



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :

Mercury Sphygmomanometer, Mobile

Product Code :

SAMCOMEDICAL6001A017

Description :

Mercury Sphygmomanometer, Mobile

Technical Specification :

Mercury Sphygmomanometer, Mobile, stand on aluminum tube with 5 plastic castors, 2 Tube adult size bladder, standard latex bulb, inflation valve, standard end valve, short latex tube with plastic connector. Mercury Sphygmomanometer, Mobile, Diagnostic Equipments Exporters, Diagnostic Equipments Manufacturers



Samco Medical

Website: www.samcomedical.com,