



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Finger Cot

Product Code :
SAMCOMEDICAL100004A002

Description :

Finger Cot

Technical Specification :

Finger Cot Size : S, M, L Feature : V shaped malleable Aluminium splints designed for easy application Poly foam padding for comfort Fits any finger, meant to keep straight for given period of time. Application : Fractures of distal inter phalangeal joints Collateral ligament injury of the distal inter phalangeal joint. Protection from Nail bed injuries of fingers



Samco Medical

Website: www.samcomedical.com,