



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Frog Splint

**Product Code :**  
SAMCOMEDICAL100004A003

**Description :**

Frog Splint

**Technical Specification :**

Frog Splint Size: S, M, L Feature: Foam padded aluminium splint in unique shape to cover the dorsal surface of finger to remain open Four extended arms to ensure good fitting. Application: Hyper flexion injuries Swan neck deformity To hold distal inter phalangeal joints in correct alignment Collateral ligament injuries of distal inter phalangeal joint



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),