



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Finger Extension Splint

Product Code :
SAMCOMEDICAL100004A005

Description :

Finger Extension Splint

Technical Specification :

Finger Extension Splint Size : Universal Feature : Malleable aluminium used to keep interphalangeal joints in right position Polyfoam padding for comfort Fits any finger Can be moulded by bare hands Application : Extensor tendon injuries of the finger Collateral ligament injuries of the inter phalangeal joint Boutonniere's deformity Immobilization to the fractured finger



Samco Medical

Website: www.samcomedical.com,