



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Wrist And Forearm Splint

Product Code :
SAMCOMEDICAL100004A009

Description :

Wrist And Forearm Splint

Technical Specification :

Wrist And Forearm Splint Size : S, M, L, XL (Left, Right) Feature : Pre shaped splint for immobilization in dorsiflexion. High quality fabric and foam padding for comfort wear. Thumb flap with velcro closure. Application : Strain and sprain in the wrist. Early cast removal. Minor fractures of the wrist and distal forearm



Samco Medical

Website: www.samcomedical.com,