



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Fracture Bracing Joint

**Product Code :**  
SAMCOMEDICAL100008A001

**Description :**

Fracture Bracing Joint Size: Universal Features: Flexible high-density unbreakable plastic. Suitable to be used on cast. Attaches directly to the above knee and below knee sections of cast. Applications: Fracture bracing to provide functional mobility and early ambulation in post acute fracture. Post acute fracture application to provide mobility and increased muscle activity. Provides early ambulation with fractures of distal femur and proximal tibia



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),