



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :

Clavicle Brace

Product Code :

SAMCOMEDICAL100010A003

**Description :**

Clavicle Brace Size : S, M, L Feature : Foam with stockinet for better comfort Straps are provided with metal fasteners. Triangular pattern for two way adjustment at back. Application : Management of clavicular injuries. Acromio-clavicular and scapular separations due to sprain Posture correction in case of stooped shoulders



Samco Medical

Website: www.samcomedical.com,