



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Samcomedical Under Arm Crutches Aluminium

**Product Code :**  
SAMCOMEDICAL100009A028



**Description :**

Samcomedical Under Arm Crutches Aluminium

**Technical Specification :**

Under Arm Crutches Aluminium Made of Aluminium. Height Adjustable. Large Size: 134-154 cm. Medium Size: 134-144 cm. Small Size: 94-114 cm



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),