



Email : sales@samcomedical.com  
Phone: 1-300-400-8211

**Product Name :**  
MUSCLE STIMULATOR

**Product Code :**  
SAMCOMEDICAL210002A001

**Description :**

MUSCLE STIMULATOR The IC based Diagnostic cum Therapeutic unit, suitable taking S.D. Curve is a marvel of electronics. The combination unit produces four different types of stimulating current i.e Galvanic, Interrupted Galvanic ,Faradic & Surged for S.D. Curve plotting and theraprutic use Output is smoothly adjustable from 0-135 volts. Specifications: GLAVANIC. : continuous well filleded D.C upto 135 volts for a maximum of 0-75mA load. INERRUPTED GALVANIC. : Pulse durations -0.01,0.03,0.1,0.3,1,3,10,30,100,300 m sec. FREQUENCY. : 6,18,30,PER MINTUE 1,3,5,10,30,50,100 per sec. SURGED-FARADIC : Faradic pulse of 0.6 m sec durations. Contraction durations: 2-4 seconds. FARADIC : Comprises of 80 pulse per second, each pulse of 0.6 ms width. ACCESSORIES : Inactive & Active Plate Electrodes.3, Ball tip Electrode with Handle.2 pair, Electrode WIRES



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),