



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
INTERFERENTIAL THERAPY UNIT

Product Code :
SAMCOMEDICAL210002A002



Description :

INTERFERENTIAL THERAPY UNIT Interferential current stimulation is effective upon the Nerves, Muscles, Arteries and Veins. The current is able to penetrate tissues to a much greater depth, with a soft stimulative effect, thus assuring an exclusive range of applications. SPECIFICATIONS: Carrier Frequency : 2 KHz and 4 KHz. Base AMF : 0 to 150 Hz. Sweep AMF : 0 to 100 Hz. Spectrum : 3 Position 1) Rectangular 1/1 s. 2) Trapezoidal 1/5 S. 3) Linear 6/6 S. Therapy Mode : 2 pole, 4 pole vector Intensity : 0-100 MA continuously adj. Intensity Balance : Available. Out-put : 3 digit indicator continuously indicates real patient current. Timer : 0 To 99 minutes digital Timer



Samco Medical

Website: www.samcomedical.com,