



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
ULTRASOUND THERAPY UNIT

Product Code :
SAMCOMEDICAL210002A003



Description :

ULTRASOUND THERAPY UNIT Micro-Processor based, computerised ultrasound therapy unit gives you a choice of treatment modes Continuous or in Three Pulsed wave forms. The compact digital unit with built in timer and friendly control panel gives a wide range of treatment situations. SPECIFICATIONS: Programme.: 24 Preset Programmes. U.sound Frequency. : 1MHz. Output Mode : Continuous & Pulsed. Variable Pulsed Frequency. : 16Hz, 48Hz, 99Hz. Intensity, Continuous Mode, Pulsed Mode : 0 to 2 watts/cm sq., 0 to 3 watts/cm sq. Treatment Head : One. Effective Area of Radiation : 5 sq. centimeters. Treatment Timer : 1 to 60 minutes, with buzzer at the end of time. Accessories. : Ultra-sound Gel bottle. Power Required : 220V.AC, 50 Hz



Samco Medical

Website: www.samcomedical.com,