



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
SHORTWAVE MEDICAL DIATHERMY-500W

**Product Code :**  
SAMCOMEDICAL210002A005



**Description :**

SHORTWAVE MEDICAL DIATHERMY-500W High power medical diathermy incorporating a high power transmitting oscillator valve & solid state rectifiers; provides R.F. power upto 500 watts for intense deep tissue heating with out discomfort at the surface. Housed in an elegant MOBILE CABINET & designed to meet heavy work load in a hospital. SPECIFICATIONS: Frequency : 27.12MHz. Output. : Continuous, maximum: 500watts. Intensity. : Adjustable, in 5 steps. Electrodes. : One pair Condenser pads with felt spacer & cover. Size: 15cm x 25cm. Tuning. : Manual, with tuning control knob. Timer. : 30 minutes timer with Auto-cut. Cooling. : With 10cm dia. instrumental fan. Voltage. : 220-224V. Accessories. : a, One, pair condenser Pads cable. b, One, Power Main cable. c, Spare Fuses. d, One, Shortwave Test Tube. e, Two Velcro Strap to tie/ hold condenser pads with the body. f, One, automatic voltage stabiliser



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),