



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
TRANSCUTANEOUS NERVE STIMULATOR

Product Code :
SAMCOMEDICAL210002A007



Description :

TRANSCUTANEOUS NERVE STIMULATOR Imported, light weight, pocket model, Two channel, Digital Tens, works on three AA battery cells. It utilizes solid state circuit to provide maximum reliability and stability.
SPECIFICATIONS: Six different stimulation modes LCD display & digital operation. Intensity adjustable to 0-9. Two channel output. Preset 10 minutes Timer. Large buttons for easy operation. Maximum output 200 Hz for severe pain



Samco Medical

Website: www.samcomedical.com,