



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
FOOT DROP BOOTS

Product Code :
SAMCOMEDICAL210003A001

Description :

FOOT DROP BOOTS A unique device to give quadriceps exercises to a person whose ankle is weak or frail. Moulded aluminium construction similar to quadriceps boot, but provided with extra support and straps at the back to be tied to the leg of the person to avoid foot drop. Weight approx. 1 kg. Comes complete with bar and collars but without weights



Samco Medical

Website: www.samcomedical.com,