



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
NIRMAL ANKLE AND LEG EXERCISER

**Product Code :**  
SAMCOMEDICAL210003A002



**Description :**

NIRMAL ANKLE AND LEG EXERCISER Aluminium moulded special designed unit, mounted on cast iron base. Unit is adjustable in height and can permit a full range of motion in eversion, inversion, dorsi flexion and plantar flexion. Resistance can be applied to these motions by addition of weights to the weight pans. Comes with a set of weights (6kg.) in disc form. Aluminium parts are bright polished and other parts finished in stove enamel



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),