



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
HEEL EXERCISER

Product Code :
SAMCOMEDICAL210003A004



Description :

HEEL EXERCISER Specially designed unit for exercising the heel & tarsus muscle groups. Unit is fitted on a laminated wooden base and patient can exercise while sitting on chair. The resistance can be increased or decreased by shifting the foot holder up or down on the lever. The pressing rod is marked and fitted with a marker for recording the improvement



Samco Medical

Website: www.samcomedical.com,