



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
ANKLE EXERCISER

Product Code :
SAMCOMEDICAL210003A005



Description :

ANKLE EXERCISER Painted aluminium boot with strap to hold the foot is mounted on a laminated base with Stainless steel sheet bracket having Two springs for ankle (Dorsi and Plantar flexion) exercises. Resistance can be increased further by tightening spring's nut. Degree of flexion can be measured and recorded



Samco Medical

Website: www.samcomedical.com,