



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
PARALLEL WALKING BAR

Product Code :
SAMCOMEDICAL210003A006



Description :

PARALLEL WALKING BAR Length 3.60 meters, with Hand rails made of iron tubing. Height and width adjusting arms made of chrome plated steel tubing of 37mm outer dia. Outer uprights made of thick steel tubing (50 mm dia). Mounted on hard-wood platform with inclination on both ends and covered with non-slip matting. Height is adjustable from 76cm to 100cm and Width from 40cm to 65cm. Adjustment can easily be done by means of setting pins at the required holes and can further be stabilised by tightening the screws Comes with detachable polished hard-wood abduction board, which fits in the holes in the middle of the platform. Steel uprights are oven baked finish



Samco Medical

Website: www.samcomedical.com,