



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
POSTURAL TRAINING MIRROR

Product Code :
SAMCOMEDICAL210003A009



Description :

POSTURAL TRAINING MIRROR Best Indian mirror, size 180cm x 60cm, protected by 5cm wide hardwood polished frame with back covered by ply board, is fitted with a pair of angular and tubular stand. Frame is mounted on 4 ball-bearing rubber castors of 50mm dia. Approximate overall size is 210cm x 75cm. Metal parts oven baked finish



Samco Medical

Website: www.samcomedical.com,