



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
EXERCISE STAIRCASE STRAIGHT TYPE

Product Code :
SAMCOMEDICAL210003A012



Description :

EXERCISE STAIRCASE STRAIGHT TYPE Sturdily built in two sections for installation in a straight line. Step arrangement : Four 6" steps on one side leading to a platform of 24" x 24" First step moves into the 2nd to make it a bus step of 12" height. Six 4" steps on the other side. Steps are 24" wide and 9" deep. Two sets of hardwood hand rails at different height to accommodate adult and children. Steps and platform made of commercial board 12 mm thickness covered with corrugated rubber matting. Sides and supporting frame of steps made of commercial board of 18 mm thickness. Hand rails and supporting bars made of teak wood, polished natural wood finish



Samco Medical

Website: www.samcomedical.com,